



2054

Limited Availability



2056



2055

Men's Size Range: from 7 to 14, half from 7.5 to 11.5

Women's Size Range: from 6 to 10, half from 5.5 to 9.5

*Mens Sizing*

|      |    |      |    |     |      |     |    |      |    |      |      |    |    |
|------|----|------|----|-----|------|-----|----|------|----|------|------|----|----|
| USA  | 7  | 7.5  | 8  | 8.5 | 9    | 9.5 | 10 | 10.5 | 11 | 11.5 | 12   | 13 | 14 |
| UK   | 6  | 6.5  | 7  | 7.5 | 8    | 8.5 | 9  | 9.5  | 10 | 10.5 | 11   | 12 | 13 |
| EURO | 40 | 40.5 | 41 | 42  | 42.5 | 43  | 44 | 44.5 | 45 | 46   | 46.5 | 48 | 49 |

*Womens Sizing*

|      |     |    |      |    |      |    |     |      |     |    |      |
|------|-----|----|------|----|------|----|-----|------|-----|----|------|
| USA  | 5.5 | 6  | 6.5  | 7  | 7.5  | 8  | 8.5 | 9    | 9.5 | 10 | 10.5 |
| UK   | 3.5 | 4  | 4.5  | 5  | 5.5  | 6  | 6.5 | 7    | 7.5 | 8  | 8.5  |
| EURO | 36  | 37 | 37.5 | 38 | 38.5 | 39 | 40  | 40.5 | 41  | 42 | 42.5 |

**VALENCIA SHOE**

2054: Valencia – Men's Mahogany

2055: Valencia – Women's Bone

2056: Valencia – Men's Bone

**MATERIALS**

Sole: Rubber

Midsole: EVA

Upper: Genuine water-resistant-treated leather and synthetic materials.

**FEATURES**

- **Harken's Radial Traction Zones:** The full 360 degree pattern of siping allows siping flexion of the sole in all directions. This pattern is placed under the 3 main pressure zones of your foot to give each area independent traction. Used in combination with a special rubber compound, traction is greatly increased.
- Integrated tongue wing to keep tongue in place.
- Neoprene lace protector.
- Lace Stopper for no-tie fastening (2056 only).
- Nonmarking, nonrusting components.
- Antibacterial silver-woven footbed kills odor-causing bacteria.
- Full-length compression-molded EVA midsole for shock absorption, neutral pronation, and all-day comfort.
- Comes with extra laces in alternative color.
- Stability Track controls ankle and foot rotation while maintaining freedom of movement.
- Metatarsal support provides all-day forefoot support in a wider toe box for improved balance.
- Active arch support reduces strain on plantar fascia and aligns knees.
- Molded heel cup provides shock absorption, support, reduces fatigue.